



**Keratinocyte Growth Supplement-animal component free  
(KGS-acf)**  
Catalog #2172

**Product Description**

Keratinocyte Growth Supplement-animal component free (KGS-acf) is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for the culture of normal keratinocytes. The supplement is formulated (quantitatively and qualitatively) to provide an optimally balanced environment that selectively supports the growth of normal keratinocytes *in vitro*.

**Product Use**

KGS-acf is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

**Storage**

Store KGS-acf at -20°C; once added to medium, store at 4°C and protect from light. Do not refreeze after thawing.

**Shipping**

Dry ice.

**Instructions for Use**

For best results, use KGS-acf in conjunction with Keratinocyte Medium-animal component free-b (KM-acf-b, Cat #2121-b).

Thaw KGS-acf at 37°C. Gently tilt the tube several times to ensure complete mixing. Spray the medium bottle and supplement tube with 70% ethanol and wipe to remove excess liquid. In a sterile field, remove the caps without touching the interior threads with fingers. Add KGS-acf to the medium and mix well. Since several components are light-labile, the medium should not be exposed to light for extended periods. We do not recommend warming medium in a 37°C water bath prior to use. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

*Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.*